

ELECTROSPINNING

- Definition:

Electrospinning is a material fabrication technique that uses an electric charge to create extremely fine fibers, on the order of nanometers. In this technique, a polymeric liquid is subjected to high electrical stress, which produces a repulsive force that stretches and stretches the liquid until a cone-shaped current is formed. When the tension is high enough, the liquid decomposes into very fine fibers that solidify upon contact with a collection surface. These fibers can be used in a wide variety of applications, such as creating materials for tissue regeneration, filter production, electronics manufacturing, and many more.

+

•

○

ELECTROSPINNING INGREDIENTS


- Recycled paper
- PVA (polyvinyl alcohol) polymer solution
- Fine needle
- 10 mls syringe
- INNTEKIE electrospinner equipment
- High voltage generator
- Collector
- Safety equipment: Safety glasses, gloves and lab coats should be used to protect themselves during the process.



ELECTROSPINNING

Time for the process

- Total process time 6 hours
- Dilution of PVA (polyvinyl alcohol) in cold 40 minutes
- Homogenization of PVA (polyvinyl alcohol) in hot 30 minutes
- Cooling of the homogenized solution 30 minutes
- Enlistment of the INNTEKIE Electrospinner equipment, syringe and receiver paper 15 minutes
- Ignition high voltage generator and power calibration 10 minutes
- Electrospinning process (nano fiber manufacturing) 4 hours



1 hour of processing



4 hour of processing





inntekie
light your ideas

INNOVATION AND TECHNOLOGY INC

ELECTROSPINNING ASSEMBLY




Universidad
Industrial de
Santander

Step by step

- 1) Weigh 5 grams of PVA.
- 2) Measure 50 mls of water
- 3) Dilute the PVA in cold water (for 40 minutes)
- 4) Heat the solution to 80°C (for 30 minutes) until the crystals are completely diluted.
- 5) Let the solution cool for 30 minutes
- 6) Take 6 milliliters with the syringe (Without using the needle)
- 7) Place the needle into the syringe and press the plunger until 2 drops of solution come out.
- 8) Put the recycled paper on the metal base of the Electrospinner INNTEKIE box
- 9) Take the syringe to the electrospinning box in the hole that best fits
- 10) Connect the negative pole to the metal plate screw
- 11) Connect the positive pole to the needle
- 12) Close the lid
- 13) Connect the high voltage generator and turn it on (red button)
- 14) Turn the knob halfway
- 15) Check the needle every hour





innTekie
light your ideas




INNOVATION AND TECHNOLOGY INC

ELECTROSPINNING ASSEMBLY



Universidad
Industrial de
Santander

PVA Dilution



Preparation of Electrospinner INNTEKIE equipment

Place metal plate
on the base

Connect pole to
ground (black
color)

Place aluminum
foil

Uncap needle


Locate syringe in
hole

Place recycled
paper on top of
aluminum foil

Connect positive
pole in the needle
(red), Purge
needle

Close lid

Turn on high
voltage
equipment and
turn knob halfway



Spinning

